


DermaSonic

More precise targeting.
Dynamic wave technology.
Visible skin-quality improvement.

Local Dynamic Micro-Massage technology to deliver four selectable frequencies — 1, 3, 4, and 10 MHz — precisely within the dermis.



MULTI-FREQUENCY TECHNOLOGY

1, 3, 4, and 10 MHz for tailored treatments.



DYNAMIC WAVE TECHNOLOGY

Delivers precise energy for enhanced results.



TARGETS MULTIPLE SKIN LAYERS

Stimulates collagen and improves skin elasticity.



VISIBLE RESULTS

Lifts, firms, and renews for smoother, healthier-looking skin.

DermaSonic



More Precise Cellular Targeting



Dynamic Wave Technology




Visible Skin-Quality Improvement



What Is DermaSonic?

Unlike single-frequency ultrasound, DermaSonic rapidly alternates between two frequencies to create pressure gradients that target the dermis, extracellular matrix, and cell groups more precisely.



What can DermaSonic Treat?


This cross-frequency energy is designed to work precisely where aesthetic results happen.

- ✓ Phonophoresis (Active Ingredient Delivery)
- ✓ Hyaluronic acid stimulation
- ✓ Supporting collagen synthesis
- ✓ Dermatitis & eczema
- ✓ Extracellular matrix balance
- ✓ Stretch marks & cellulite



The Result? A non-invasive treatment that lifts, firms, and renews skin from within — with no downtime.

Treatment Probes & Frequencies



Probe	Skin Layer	Frequency
 Facial Probe	Dermis Surface	3 / 10 MHz
 Body Probe	Body / Deep Skin	1 / 3 MHz
 Eye Probe	Mid-Shallow Dermis (Eye)	4 MHz (RF)
 Cold & Hot Hammer	All Layers + Surface	-15 / +55°C